

Supplementary **Table 1.** Primers used in real-time PCR

Gene	Sense (5'-3')	Antisense (5'-3')	Product Length
<i>JAG1</i>	GCGTCATTGTGTTACCTGCG	GAGCCCGGCCTCCTTTTATT	227
<i>NOTCH1</i>	TGAATGGCGGGAAGTGTGAA	ATAGTCTGCCACGCCTCTG	161
<i>VEGF</i>	GAGGAGCAGTTACGGTCTGTG	TCCTTTCCTTAGCTGACACTTGT	96
<i>MMP9</i>	TGTACCGCTATGGTTACACTCG	GGCAGGGACAGTTGCTTCT	97
<i>GAPDH</i>	ACATCGCTCAGACACCATG	TGTAGTTGAGGTCAATGAAGGG	143
<i>ACTB</i>	TGCCCATCTACGAGGGGTATG	TCTCCTTAATGTCACGCACGATTT	156